

# Defibrillationselektroden

## Klebeposition



### Reanimation

Anterior/Posterior

und

### Stimulation

(Pacing)





Elektroden:

8900-0214

**OneStep CPR Complete** Elektrode

### Reanimation

Anterior/Anterior



Elektroden:

8900-0214

**OneStep CPR Complete** Elektrode

Oder



8900-0217

**OneStep CPR A/A** Elektrode

### Stimulation

(Pacing)

3 Lead EKG nötig!



Elektroden:

8900-4003

**Stat Padz** Elektrode



oder

8900-4005

**Pro Padz** Elektrode (X-Ray)

### Kardioversion

Anterior/Posterior



Elektroden:

8900-2303

**Pro Padz (Liquid)** Elektrode

oder

8900-2105

**Pro Padz (X-Ray, Liquid)** Elektrode